



Tips for Successful HFA Inhaler Use

HFA Inhalers	Contact Information	Priming	Special Instructions	Cleaning
		<ul style="list-style-type: none"> - All priming or spraying should be done away from face. - Shake well and then: 		<ul style="list-style-type: none"> - Remove canister and cap from mouthpiece. - Never immerse canister in water!
Pro Air™ HFA (albuterol sulfate)	Teva Specialty Pharmaceuticals LLC www.proairhfa.com 1.888.482.9522	- Spray 3 times when new or when not used for more than 2 weeks	None	<ul style="list-style-type: none"> - Wash mouthpiece with warm running water for 30 seconds <u>once a week</u>. - Shake to remove excess water. Air dry overnight.
Proventil® HFA (albuterol sulfate)	Schering Corporation www.proventil.com	- Spray 4 times when new or when not used for more than 2 weeks	None	<ul style="list-style-type: none"> - Wash mouthpiece with warm running water for 30 seconds <u>once a week</u>. - Shake to remove excess water. Air dry overnight.
Ventolin® HFA (albuterol sulfate)	GlaxoSmithKline www.ventolin.com 1.888.825.5249	<ul style="list-style-type: none"> - Spray 4 times - when new or when not used for more than 2 weeks, or - <u>when inhaler has been dropped</u> - Shake well between sprays. 	<ul style="list-style-type: none"> - Store at room temperature (68°F to 77°F) - Ventolin HFA has a counter to track doses. Discard inhaler when counter reads 000. - Store with mouthpiece down. 	<ul style="list-style-type: none"> - Wash mouthpiece with warm running water for 30 seconds <u>once a week</u>. - Shake to remove excess water. Air dry overnight. - Prime once after cleaning.
Xopenex® HFA (levalbuteral tartrate)	Sepracor Inc. www.xopenex.com 1.888.394.7377	- Spray 4 times when new and if not used for 3 days.	<ul style="list-style-type: none"> - Store at room temperature. - Prevent direct sun light. - Store with mouthpiece down. 	<ul style="list-style-type: none"> - Wash mouthpiece with warm running water for 30 seconds <u>once a week</u>. - Shake to remove excess water. Air dry overnight.

Each inhaler contains 200 sprays